


Tail Lysis Buffer

Tail Lysis Buffer is ready to use solution that enables simple genotyping procedure.


- » Ready-to-use solution
- » DNase, RNase free

Application 1: Genotyping of mouse tail


Procedure



Result 1



Result 2



Sry gene (132 bp), which only exists in male mouse Y-chromosomes was increased in male derived DNA, but was not increased in female derived DNA. The result shows that PCR operates efficiently.

Ordering Information

Product Name	Storage	Product No.	PKG Size
Tail Lysis Buffer	RT	06169-95	500 ml

Related Products

Product Name	Storage	Product No.	PKG Size
Proteinase K Solution (>700U/ml)	R	04130-06	2 ml
		04130-64	10 ml
Phenol:Chloroform:Isoamyl Alcohol 25:24:1 Mixed, pH7.9	R	25970-14	100 ml
		25970-56	400 ml
2-Propanol	RT	03065-35	500 ml
Ethanol (99.5)	RT	08948-25	500 ml

[Storage] RT = Room temperature, R = Refrigerator