


Silver Staining Kit

Sliver staining method is high sensitive method for detecting nucleic acids and proteins in polyacrylamide gel. We offer three types of silver staining kits, each having unique features for your experimental needs.

Comparison of Each Property

	Sil-Best Stain One	Sil-Best Stain-Neo	Sil-Best Stain
2-Dimensions	Excellent	Poor	Good
SDS-PAGE	Good	Good	Good
Nucleic Acid	Poor	Good	Fair
Step	12	6	14
Staining time	80 min.	60 min.	110 min.

Comparison of Each Procudre



Comparison of Each Staining Image

